Nonisothermal Kinetics Analysis of the Dehydration of Ziprasidone Hydrochloride Monohydrate by Thermogravimetry

A. RAVIKIRAN*, M. ARTHANAREESWARI, P. KAMARAJ, AND C. PRAVEEN Department of Chemistry, SRM University, SRM Nagar, Kattankulathur-603 203, Kancheepuram, India

Ravikiran, et al.: Dehydration Kinetics of Ziprasidone Hydrochloride Monohydrate

In the current work the kinetics of dehydration of ziprasidone hydrochloride monohydrate was studied by nonisothermal thermogravimetry. Ziprasidone hydrochloride monohydrate was heated from 30 to 150° with a heating rate of 5° per min under nitrogen gas atmosphere and weight loss data were collected. Powder X-ray difraction was used to characterize the solid before and after dehydration. The well accepted Coats-Redfern model fitting approach was applied to the thermogravimetry data for the kinetic analysis. Thirteen solid state reaction models were studied; among them one-dimensional diffusion model was found to be the best fit model for this reaction with an excellent correlation 0.9994. The Arrhenius parameters, activation energy, and pre-exponential factor were determined, the values were found to be 28 k.cal/mol and 9.53×10^{13} sec⁻¹, respectively.

Key words: Model fitting, kinetics, dehydration, thermogravimetry, ziprasidone hydrochloride

Ziprasidone hydrochloride (ZH), chemically 5-[2-[4-(1,2benzisothiazol-3-yl)-1-piperazinyl] ethyl]-6-chloro-1,3-dihydro-2H-indol-2-one, monohydrochloride, represented in fig. 1, is a typical antipsychotic drug used for the treatment of schizophrenia, a mental disorder. It is commercially available as Geodon[®] capsules^[1]. Anhydrous and monohydrate forms of ZH are known from literature^[2].

In view of increasing regulatory concern related to solid form of the drug substance in the recent days, it is highly important to understand the physical stability of different solid forms of a drug substance and their inter-convertibility^[3,4]. Understanding the kinetics of solid state reactions in relation to

phase transformations of drug substance leaves an opportunity to optimize the operational conditions so as to prevent and have control over the changes in solid from, during manufacturing process and shelf life of the drug product. In the present work dehydration of ZH monohydrate to ZH anhydrous form, as represented in Eqn. (1) was studied by thermogravimetry (TG). ZH.H₂O \rightarrow ZH+H₂O.(1)



*Address for correspondence E-mail: ravianalytical@gmail.com

Fig. 1: Molecular structure of ziprasidone hydrochloride monohydrate.

TG is often used to study the kinetics of solid state reactions involving weight change during reaction. Some examples of such reactions include thermal decomposition, oxidation, and dehydration or desolvation^[5,6]. Since the number of data points is high enough for kinetic analysis of weight loss data, TG has become useful for this purpose. In the present work TA's Thermo gravimetric analyzer (TGA), Model Q5000IR was used for the study. Temperature calibration of TGA was performed by measuring Nickel Curie point (C_p) at 354°. The Curie point is a physical property of a metal; it is the temperature point above which the metal loses ferromagnetic properties. Weight calibration of TGA was performed by using 50 mg standard reference material. About 15.7 mg of ZH-monohydrate was taken into aluminum pan. The instrument was programmed to heat the specimen from 30 to 150° with a heating rate 5° per min. Nitrogen was used as a purge gas with a flow rate of 50 ml/min. Model fitting approach^[7-9] was applied on the collected weight loss data using Coats-Redfern method^[10]. Powder X-Ray diffraction (PXRD) patterns of the ZH-Anhydrous and ZH-Monohydrate were collected by using Bruker D8 ADVANCE Powder X-ray Diffractometer equipped with Cu anode.

General form of rate equation of solid state reactions used for isothermal analysis is $\frac{d\alpha}{dt} = k(T)f(\alpha) \dots (2)$, where α is the extent of reaction, t is the time, T is the temperature, and $f(\alpha)$ is reaction model. For the current dehydration study: $\alpha = \frac{w_o - w_t}{B.w_o} \dots (3)$, where w_o , w_t and B are the initial mass of the sample, the mass at time t and the fraction of weight loss for complete dehydration of the reaction, respectively. Integral form of Eqn. (2) is $g(\alpha) = kt$, with $g(\alpha) = \int_0^\alpha \frac{1}{f(\alpha)} d\alpha \dots (4)$. According to Arrhenius equation $k = Ae^{\frac{-Ea}{RT}}$ $\dots (5)$, where k, E_a , and A are the rate constant, the activation energy, and the pre-exponential factor, respectively. $g(\alpha) = Ae^{\frac{-Ea}{RT}}t \dots (6)$, under nonisothermal conditions integral form of Eqn. (6) is $g(\alpha) = \frac{A}{\beta} \int_0^T e^{\frac{-Ea}{RT}} \dots (7)$, where β is the heating rate. If $\frac{Ea}{RT}$ is replaced by a variable x, the temperature integral Eqn. (7) becomes $g(\alpha) = \frac{AEa}{\beta R} \int_x^{\infty} \frac{e^{-x}}{x^2} dx$

....(8). If
$$p(x) = \int_{x}^{\infty} \frac{e^{-x}}{x^{2}}$$
, then the Eqn. (8) becomes $g(\alpha) = \frac{AEa}{\beta R} p(x) \dots (9)$, where, $p(x)$ is the

exponential integral, it has no analytical solutions but has many approximations^[11]. The Eqn. (9) can be used for several model fitting methods for the analysis of nonisothermal kinetics.

Model-fitting approach involves fitting different models to α -temperature curves. A model is a mathematical expression, developed based on mechanistic assumptions, transforms a solid state reaction process into a rate equation. Therefore, different rate expressions are produced from these models. To explain solid state reactions, there are four models such as (i) nucleation, (ii) geometrical contraction, (iii) diffusion, and (iv) reaction order models. Briefly, nucleation model is often used for reactions like crystallization, crystallographic transition, decomposition, adsorption hydration, and desolvation. The geometrical contraction model assumes that nucleation occurs rapidly on the surface of the crystal and the rate of degradation is controlled by the resulting reaction interface progress toward the center of the crystal. The diffusion model best explains solid state reactions involving gaseous products, where the reaction rate is controlled by the movement of reactants or products from the reaction interface or product layer. The order-based models are the simplest models as they are similar to those used in homogeneous kinetics. In these models, the reaction rate is proportional to concentration, amount or fraction remaining of reactant(s) raised to a particular power (integral or fractional) which is the reaction order^[12].

In the current study popular Coats-Redfern method was used in which the asymptotic series expansion for approximating Eqn. (9) is used to get the following equation: $\ln \frac{g(\alpha)}{2} = \ln \left[\frac{AR}{2RT}\left(1-\frac{2RT}{2RT}\right)\right] - \frac{Ea}{2RT}$(10),

where
$$\overline{T}$$
 is the mean experimental temperature. Using

this Eqn. (10), the values of E_a and A can be obtained from slope and intercept values, respectively, from the graph plotted for $ln \frac{g(\alpha)}{T^2}$ versus 1/T for different models, as represented in Table 1.

The TG data indicated a total weight loss of about 4.067% w/w in the temperature range from about

DEHYDRATION OF ZIPRASIDONE HYDROCHLORIDE MONOHYDRATE USING COATS-REDFERN METHOD					
Reaction model	$f(\alpha)$	g(a)	Ea (kcal/mol)	A	-R
Nucleation models					
Power law	4 α ^{3/4}	$\alpha^{1/4}$	2	0.01	0.9988
Power law	3 α ^{2/3}	α ^{1/3}	4	0.03	0.9990
Power law	2 α ^{1/2}	α ^{1/2}	0	0.00	0.9945
Avrami-Erofeev	4 (1-α) [-ln (1-α)] ^{3/4}	[-ln (1-α)] ^{1/4}	3	0.01	0.9956
Avrami-Erofeev	3 (1-α) [-ln (1-α)] ^{2/3}	[-ln (1-α)] ^{1/3}	4	0.12	0.9963
Avrami-Erofeev	2 (1-α) [-ln (1-α)] ^{1/2}	[-ln (1-α)] ^{1/2}	7	9.18	0.9969
Diffusion models					
One dimensional diffusion	1/2 α ⁻¹	α^2	28	9.53×1013	0.9994
Diffusion control (Janders)	2 (1-α) ^{2/3} [1-(1-α) ^{1/3}] ⁻¹	$[1-(1-\alpha)^{1/3}]^2$	31	1.47×10 ¹⁵	0.9983
Diffusion control (Crank)	$3/2[(1-\alpha)^{-1/3}-1]^{-1}$	1-2/3 α -(1-α) ^{2/3}	30	2.66×10 ¹⁴	0.9987
Reaction order and geometrica	l contraction models				
Mampel (first order)	1-α	[-ln (1-α)]	16	3.06×10 ⁶	0.9974
Second Order	(1-α) ²	(1-α) ⁻¹ -1	18	2.20×10 ⁸	0.9944
Contracting cylinder	2 (1-α) ^{1/2}	1-(1- α) ^{1/2}	14	2.12×10 ⁵	0.9985
Contracting Sphere	3 (1-α) ^{2/3}	1-(1-α) ^{1/3}	15	2.70×10⁵	0.9981

 TABLE 1: SOLID STATE REACTION MODELS USED AND ARRHENIUS PARAMETERS FOR NONISOTHERMAL

 DEHYDRATION OF ZIPRASIDONE HYDROCHLORIDE MONOHYDRATE USING COATS-REDFERN METHOD

One-dimensional diffusion model has resulted best correlation (0.9994); The calculated values of the activation energy (Ea) and pre-exponential factor (A) of this model were considered for the current solid state reaction



Fig. 2: Thermogravimetric analysis graph of ZH-monohydrate. The curve represents significant weight loss in the temperature range from about 50° to 80°, about 50% reaction completed at about 73°.

30 to 90° and it was also observed that the fraction of extent of reaction was reached 0.5 at about 73°. The TG thermogram is represented in fig. 2. The comparison of PXRD patterns of ZH-anhydrous and ZH-monohydrate is represented in fig. 3, which shows the change in the PXRD profile of ZH-monohydrate after dehydration process, indicating conversion of the monohydrate into an anhydrous form, having a different crystalline structure.

By using various expressions for $g(\alpha)$ in Eqn. (10), Arrhenius parameters were calculated from the plot of $ln \frac{g(\alpha)}{T^2}$ against *1/T*. The set of calculated Arrhenius parameters for dehydration of ZH monohydrate are



Fig. 3: Powder X-ray diffraction overlay of ZH-anhydrous and ZH-monohydrate.

ZH anhydrous (-) and ZH monohydrate (-) solid forms show distinct powder X-Ray Diffraction profiles.

listed in Table 1. For each model, goodness fit was determined by the correlation value, as mentioned in Table 1. All models have resultant correlation greater than 0.99. The first six reaction models have very small values of pre-exponential factor demonstrating that these models cannot explain the reaction kinetics. Among the remaining models, one-dimensional diffusion model has got maximum correlation of 0.9994.

In conclusion, dehydration of ZH-monohydrate was successfully studied using TG and PXRD. A distinct Powder XRD pattern was observed for dehydrated solid when compared with the monohydrate, indicating transformation into different crystalline phase. The nonisothermal weight loss data obtained from TG were used for kinetics analysis, one dimensional diffusion model was found to be the best fit among 13 models evaluated using Coats-Redfern method and the calculated values of *Ea* and *A* are 28 kcal/mol and 9.53×10^{13} sec⁻¹, respectively. Thus, reliable evaluation of kinetic parameters of dehydration process of pharmaceutical hydrates is very important to optimize the operating conditions required during manufacturing process and shelf-life of its formulation. TG in combination with powder XRD can be used to evaluate dehydration kinetics and to determine Arrhenius parameters, in relatively short time.

REFERENCES

- Geodon.com [Internet]. New York: Pfizer Inc.; c2012-Available from: https://www.geodon.com/about-geodon.aspx [Last cited on 2012 Nov 20].
- United States Patent Number 5,312,925 (1994 May 17), Monohydrate of 5-(2-(4-(1,2-benzisothiazol-3-yl)-1-piperazinyl)-ethyl)-6-chloro-1,3-dih ydro-2H-indol-2-one-hydrochloride. Available from: http://patft.uspto. gov/[Last updated on 2010 Aug 26; cited on 2012 Nov 20].
- 3. Yu LX, Furness MS, Raw A, Outlaw KP, Nashed NE, Ramos E, Miller SP,

et al. Scientific considerations of pharmaceutical solid polymorphism in abbreviated new drug application. Pharm Res 2003;20:531-6.

- Raw AS, Furness MS, Gill DS, Adams RC, Holcombe FO Jr, Yu LX. Regulatory considerations of pharmaceutical solid polymorphism in Abbreviated New Drug Applications (ANDAs). Adv Drug Deliv Rev 2004;56:397-414.
- 5. Giron D. Applications of thermal analysis in the pharmaceutical industry. J Pharm Biomed Anal 1986;4:755-70.
- MacCallum JR, Tanner J. The kinetics of thermogravimetry. Eur Polym J 1970;6:1033-9.
- Doyle CD. Kinetic analysis of thermogravimetric data. J Appl Polym Sci 1961;5:285-92.
- Vyazovkin S, Wight CA. Model-free and model-fitting approaches to kinetic analysis of isothermal and nonisothermal data. Thermochim Acta 1999;340:53-68.
- 9. Tanaka H. Thermal analysis and kinetics of solid state reactions. Thermochim Acta 1995;267:29-44.
- Coats AW, Redfern JP. Kinetic parameters from thermogravimetric data. Nature 1964;201:68-9.
- Doyle CD. Series Approximations to the equation of thermogravimetric data. Nature 1965;207:290-1.
- 12. Khawam A, Flanagan DR. Solid-state kinetic models: Basics and mathematical fundamentals. J Phys Chem B 2006;110:17315-28.

Accepted 12 April 2013 Revised 07 April 2013 Received 16 January 2013 Indian J Pharm Sci 2013;75(3):361-364